

LunaLegFlat

Guide

Loop in strap



1) Pull strap through a deflection loop with the fleece side facing upwards

2) Fasten the end of the strap on the fleece.

Attach LunaLeg on the leg



1) Place the LunaLeg on your thigh. The strap with the fleece and Velcro side facing outwards.



2) Pull the strap through the loop. Fasten the end of the strap on the fleece.

Mount Smartphone

First push the smartphone into the upper silicone band, then into the lower band. The silicone bands must be stretched with some force. Stretch the bands over all 4 corners.



This way, the smartphone is correctly and firmly in the holder:



TIP: First mount the smartphone, then attach the LunaLegFlat on your Leg.